

- RNA isolation phase (w/ [RNEasy Plus Mini Kit](#)):
 - Every step at room temperature
 - Gather sample (cell, organoid) ($< 1 \times 10^7$ cells/sample).
 - Lyse directly or trypsinized and collected as a cell pellet by adding 400 μ L Buffer RLT Plus and keep it in 1.5 mL tube
If organoid, use P1000 pipette tip to break up organoid and do a swirling motion “beating an egg.”
 - Homogenize the lysate by vortexing for 1 min
>> Pause point: -20°C or -80°C for several months
 - Transfer lysate to purple gDNA Eliminator Spin Column. Pipette up and down until foamy.
 - Centrifuge for 30s at 8000g. Save Flow-Through (/FT). Discard column.
*All centrifuge at $20-25^{\circ}\text{C}$ (not cool below 20°C)
 - Add 350 μ L of 70% EtOH to FT. Mix and immediately transfer 700 μ L of FT to the pink RNEasy Spin Column.
 - Close lid & Centrifuge for 15s at 8000g. Discard FT. Save column.
 - Add 700 μ L Buffer RW1 to the same pink column.
 - Close lid & Centrifuge for 15s at 8000g. Discard FT. Save column.
 - Add 500 μ L Buffer RPE to the same pink column.
 - Close lid & Centrifuge for 15s at 8000g. Discard FT. Save column.
 - Add 500 μ L Buffer RPE to the same pink column.
 - Close lid & Centrifuge for 2min at 8000g. Discard FT. Save column.
 - Dry column for 2min.
 - Place pink column in a new 1.5-mL Eppendorf tube for elution.
 - Add 35-60 μ L dw directly onto the membrane
 - Close lid & Centrifuge for 1min at 8000g to elute RNA.
 - Measure [RNA] with NanoDrop.
 - Store RNA at 4°C until further use.
- RT phase:
 - Dilute RNA with dw in PCR-tube strip.
 - Make RT Master Mix.
 - Distribute RT MM into diluted RNA in PCR-tube strip.
 - On PCR machine, incubate at 42°C for 2 hours then 4°C .
 - Store cDNA RT product at 4°C until further use.
 - Formula for RNA dilution per sample:
 - Ideally, target mass of RNA = 1000 ng (max 2000 ng)
 - If smallest [RNA] = 65 ng/ μ L, then set target mass at 65 ng/ μ L \times 15 μ L = 1000 ng
 - Calculate volume of RNA needed to reach target mass for each sample.
 - Top off with dw to reach total volume of 15 μ L
 - Total volume for RNA diluted = 15 μ L
 - Formula for RT Master Mix in 2-mL Eppendorf tube:
 - Multiply each volume below by (the number of reactions needed + 10% for pipetting error)

- 2 uL 10X RT Buffer (Applied Biosystem)
- 2 uL 10X RT Random Primers
- 1 uL RT Enzyme (MultiScribe Reverse Transcriptase)
- 0.5 uL dNTP
- Total volume for RT MM per sample = 5.5 uL
- Total volume for RT product per sample = 20.5 uL
- Thermocycler setup for reverse transcriptase (two-step)

	Step 1	Step 2	Step 3	Step 4
Temperature (°C)	25	37	85	4
Time	10 min	120 min	5 min	∞

- 5RT with BioRad's iScript 5x MM
 - Dilute RNA to 1000ng or 1500 ng in PCR strip → Total Volume 8uL
 - Add 2uL 5x MM to make it to 10uL
- qPCR phase (ddCT):
 - Formula for each well:
 - 5 uL 2X PowerSYBR Green MM
 - 0.2 uL 10uM primer mix
 - 0.2 uL RT product/template
 - 4.6 uL dw
 - Step 1: Plate Planning
 - Suppose: 5 samples, 10 genes
 - Mixtures: SYBR MM + Template; dw + primer (not waste more SYBR!)
 - Prepare 5 tubes in strip for 5 templates (dist. horizontally)
 - Each tube has enough for 10 → 15 rxns
 - Prepare 10 tubes in strip for 10 genes (dist. vertically)
 - Each tube has enough for 5 → 10 rxns
 - Step 2: Volume Planning
 - SYBR MM + Template Mixture:
 - SYBR MM: 5 uL x 15 rxns = 75 uL
 - Template: 0.2 uL x 15 rxns = 3 uL
 - Total = 78 uL per PCR tube
 - Distributing horizontally 5.2 uL per well
 - Dw + Primer Mixture
 - Primers: 0.2 uL x 10 rxns = 2 uL
 - Dw: 4.6 uL x 10 rxns = 46 uL
 - Total = 48 uL per PCR tube
 - Distributing vertically 4.8 uL per well
 - Step 3: Doing it!
 - Distribute with the electronic multi-channel pipette
 - Each tip has a max volume of 30 uL

- qPCR phase (Nilay 2021.8.24) :

- Make a standard

- Tube A = Pool 3 uL cDNA from each condition to make a most concentrated standard
e.g. 6 conditions -> $3 \times 6 = 18$ uL
 - Tube B = Dilute tube A in dw to make enough standard volume for each gene
e.g. 10 genes -> dilute to the final volume at least 10×4.5 uL = 45 uL
 - Serial dilution of Tube B in 1:5 manner 4 times to make tube C, D, E

	Relative cDNA
Tube B	1:1 (125x)
Tube C	1:5(25x)
Tube D	1:25(5x)
Tube E	1:125(1x)

- Formula for each rxn (10 uL/rxn)



- 5 uL 2X PowerSYBR Green MM
 - 0.2 uL 10uM primer mix (each)
 - 2 uL 10-fold dilution RT product/template
 - 2.8 uL dw

- Plate planning

- Each plate has to have 1) house keeping gene 2) standard reaction for each gene

- Thermocycler (Bio-Rad CFX)

Activation of enzyme	95c	10 min	Activate AmpliTaq Gold DNA polymerase
Denature	95c	15 sec	40x
Annealing&Elongation	60c	60 sec	
Plate readout			
Melt curve analysis	55c	0.5c	
	90c	increment every 5 sec	



[illegible][illegible][illegible]