

- **Workflow:**
 - **Generating Competent Bacteria**
 - **Digestion & Ligation (place desired insert into plasmid vector)**
 - **Bacterial Transformation**
 - **Plate Growth: Select for bacteria that uptake the vector**
 - **Colony PCR (low-resolution genotyping)**
 - **Sanger Sequencing (high-resolution genotyping)**
 - **Overnight culture (of the inoculated, genotyped, and correct clones)**
 - **Or Bacterial Freezer Stock Generation**
 - **Plasmid Purification**
 - **Repeat**

- **Generating Competent Bacteria by CaCl₂ method**
 - Strain: Invitrogen E.coli Stbl3
 - 100-200x dilution from frozen stock
 - Grow 2-3 hrs in exponential phase
 - Incubate on ice (to synchronize all cells, not mitosis phase)
 - Wash with pre-chilled 100 mM CaCl₂, while on ice
 - Store competent bacteria at -80 C
 - Expected Yield: 5e6 to 2e7 transformed colonies/ ug of plasmid DNA
 - Stage 1: Preparation of Cells
 - Pick a single colony from a plate that has been incubated overnight at 37 C.
 - Inoculate that colony into 100-mL LB broth in a 1-L flask
 - Incubate at 37 C for 3 hours
 - Begin harvesting when OD 600 reaches 0.35
 - Transfer bacteria to 50-mL conical tubes on ice for 10 min
 - Centrifuge at 2700 g for 10 min at 4 C.
 - Decant supernatant. Invert tubes on a paper towel to dry out all supernatant
 - Resuspend with 30 mL of ice-cold MgCl₂+CaCl₂ solution (80 mM MgCl₂, 20 mM CaCl₂). Vortex.
 - Centrifuge at 2700 g for 10 min at 4 C.
 - Decant supernatant. Invert tubes on a paper towel to dry out all supernatant
 - Resuspend with 4 mL of ice-cold 0.1 M CaCl₂ (2 mL per 50 mL of original culture). Vortex.
 - Dispense the competent bacteria into aliquots and store at -80C

- **Generation of Insert** (Oligo can be encoding for gRNA to be inserted into H3 or pXPR_003 plasmid vectors)
 - Go to UCSC Genome Browser
 - Type in Gene Name that you want to Knock-Out
 - Scroll down to “Genes and Gene Prediction” section

- Display “CRISPR Targets” in “full”
 - The green sequences are better NGG
 - Order oligo synthesis from IDT
 - Find 20210621_gRNA_design.xlsx to use as template
 - Add "g" to 5' end of TOP strand
 - Complete TOP strand: cacc+g+seq
 - Reverse Complement of TOP to get BOTTOM
 - Add “aaac” to 5' end of BOTTOM strand
 - Name the ssDNA oligo to be purchased
 - Liquidify oligos
 - 100 uM per oligo: 10 nmol mixed in 100 uL dw
 - Formula for oligos mixture in PCR tube strip:
 - 5 uL TOP strand
 - 5 uL BOTTOM strand
 - Phosphorylation of 5' end of oligo (synthetic oligos do not have 5' phosphate)
 - Formula in PCR tube strip:
 - 1 uL T4 DNA ligase 10X Buffer
 - 0.5 uL PNK Enzyme (polynucleotide kinase)
 - 6.5 uL dw
 - (2) uL oligos mixture
 - Put in PCR machine:
 - 37 C for 30 min
 - 95 C for 5 min
 - 28 Cycles cooldown
 - -2.5 C decrement per cycle
 - 30 sec per cycle
 - Before ligation, dilute the phosphorylated oligos 10X: 90 uL dw + 10 uL oligos mixture
- **Digestion & Ligation:** Place desired insert into vector
 - **Digestion of vector (linearization)**
 - **Ligation**
 - Formula in PCR tube strip:
 - 1 uL T4 DNA ligase
 - 1 uL T4 DNA ligase 10X Buffer
 - 1-2 uL linearized vector H3 (enough to get 50 ng of vector per rxn)
 - 2 uL phosphorylated oligos mixture
 - Top off with dw to get final volume of 10 uL per rxn
 - Mix by pipetting up and down or tap tap
 - Incubation
 - Sticky ends: room temperature for 10 min (or 16 C overnight)
 - Blunt ends: room temperature for 2 hours (or 16 C overnight)
 - Heat inactivate at 65 C for 10 min

- Chill on ice
- **Bacterial Transformation by Heat-shock Method**
 - Get aliquots of competent E coli Stbl3 bacteria from -80 C. Thaw on ice for 15 min.
 - Distribute in Each PCR tube: 50 uL bacteria + 5 uL ligation product (plasmid DNA with insert) (10% ratio of DNA:bacteria)
 - Vortex to mix. Slight spin down with bench-top microfuge.
 - Incubate mixture on ice for 30 min.
 - Heat shock @ 42 C for 45 sec.
 - Incubate on ice for 2 min.
 - In 2-mL Eppendorf tubes, add 400 uL of pre-warmed SOC media per tube.
 - Transfer transformed bacteria into these SOC-filled Eppen tubes.
 - Keep in warm room on shaker for 2 hours for regeneration
 - Using a P1000 tip, suck up the bacteria and spread on to Cb+Lb petri dishes.
 - Leave dishes in warm room overnight.
 - Notes:
 - SOC media bought from Boston Bioproducts
 - Recipe
- **Colony PCR (low-resolution genotyping)**
 - Add 50-75 uL Lb+Cb broth per well to 96-well PCR plate (leave out the edges)
 - Pick 3 colonies per agar+Cb dish (triplicates) with sterile P20 tips
 - These colonies have picked up the plasmid
 - Incubate 1 hour in warm room 37 C
 - Formula in 384-well plate:
 - 12.5 uL ThermoScientific PCR 2X MM
 - 11 uL dw
 - 1 uL of vector-specific Primer Mix (10 uM each from 80dw+10F+10R mixture)
*OV267-F,OV268-R for pLentiCRISPRV2
 - DNA template (heat lysis of bacteria)
 - Total Volume: 25 uL
 - Keep the 96-well PCR plate, if PCR genotype is fine → use it as stock
 - Touchdown PCR Phire 1kb TD65 check machine
 - 1st 10 cycles: 98-65-72
 - 2nd 10 cycles: 98-60-72
 - Non-permissive 65 C annealing temperature → permissive 60 C annealing temperature
 - Initially, only very specific hybridizations are allowed to occur
 - Temperature increment: -0.5 C per cycle
 - Primer design:
 - Spanning U6 promoter to after the gRNA scaffold (370 bp region)

- **Gel Electrophoresis (to verify we have product)**
 - See agarose gel making protocol
- **Sanger Sequencing (high-resolution genotyping)**
 - Send pre-mixed template + primer
 - Diluted primer Formula
 - 5 uL FORWARD primer
 - 255 uL dw
 - Total Volume 300 uL
 - Formula of sample to submit to Genewiz
 - 13 uL diluted F primer
 - 2 uL PCR product
 - GeneWiz Instructions
 - PCR Product (Purified)
 - Service Type: Premix (Primer + Template)
 - "Purification": "Enzymatic"
 - Put in the bag a note: the barcode, # of samples submitted, Sanger Sequencing
- **Overnight culture (of the inoculated, genotyped, and correct clones)**
 - Prepare 50 mL Lb + Cb in 1 50-mL conical tube for every correct clone.
 - Add ~75 uL (all) of the regenerated clone from the 96-well PCR plate to each conical tube
 - Vortex.
 - Split into 2 50-mL conical tubes (duplicates)
 - Overnight 37 C shaker for 14-16 hours (no more or less)
- **Stock Generation by Freezing**
 - From overnight culture of the **inoculated, genotyped, and Sanger-correct clones**
 - Add 1m L overnight growth to 2-mL cryovial
 - 1.5-mL Eppendorf tubes can snap open at -80
 - Top off with 1 mL 50% glycerol (dilute 1:1 beforehand with dw)
 - Store at -80 C
- **Plasmid Purification by QIAGEN PlasmidPrep Midi Plus Kit; high-yield protocol**

Table 3. Maximum recommended culture volumes

Copy number	Standard protocol	High-yield protocol
High-copy plasmid*	20–25 ml	25–35 ml
Low-copy plasmid*†	50 ml	–

* For high-copy plasmids, expected yields are ▲100–200 µg for the QIAGEN Plasmid *Plus* Midi Kit using the standard protocol and ●150–250 µg for the QIAGEN Plasmid *Plus* Midi using the high-yield protocol. For low-copy plasmids, expected yields are ▲30–100 µg for the QIAGEN Plasmid *Plus* Midi Kit using these culture values with the standard protocol.

- Collect bacteria from overnight culture, suspension total volume 50 mL

- Centrifuge 50-mL Falcon tubes for 10min at 3000g, discard supernatant
- Vortex Falcon tubes to break pellets.
- Add 4 mL P1 Lysis Buffer (TEG, hypo-osmotic lysis)
 - Shake, vortex or pipet until no cell clump. Wait 3-4 min
- Add 4 mL P2 Buffer (blue, SDS detergent, basic pH) to denature bacterial & plasmid DNA
 - Carefully flip tubes 10x to precipitate bacterial gDNA, protein, lipid
 - Wait 3-4 min
 - viscous lysate, mix until homogenize; the solution should be blue
- Add 4 mL S3 binding buffer (acetic acid) to quench
 - Carefully flip tubes 10x to mix till the blue color disappear
 - white precipitate = gDNA, proteins, cell debris, SDS
 - Do not wait
- Centrifuge for 5 min at 4500g, Prepare the vacuum manifold and QIAGEN spin columns
- Insert the plunger into the QIA filter cartridge and filter the cell lysate into a new tube
 - *filter into the extension tube directly above the spin column
- Add 2 mL BB buffer to the cleared lysate, invert 4-6 times
 - *or add BB to the extension tube
- Transfer lysate to the spin column
- Switch on the vacuum source ~300 mbar
- Wash DNA with 0.7 mL ETR buffer, switch on vacuum then switch off
- Wash DNA with 0.7 mL PE buffer, switch on vacuum then switch off
- Put the column inside a 1.5 mL tube, centrifuge 10,000 g x 1 min to remove residual wash buffer
- Place a column into a new 1.5 mL tube; elute DNA with 200 uL EB buffer or water to the center of a spin column, wait 1 mins and centrifuge 10,000g x 1 min
 - *Elute with 100 uL if high concentration needed
 - *DNA stored in water may degrade over time
- Store DNA at -20*c